



























FEMORAL IMPLANT SPECIFICATIONS

CATALOG #	SIZE Ø x LENGTH [MM]	PROXIMAL FIXATION [MM]	DISTAL FIXATION [MM]
<i>Standard femoral implants (long threaded distal fixation)</i>			
FD-032(L)-SS	Ø3.2 x 202		
FD-040(L)-SS	Ø4.0 x 338		
FD-048(L)-SS	Ø4.8 x 410		
FD-056(L)-SS	Ø5.6 x 410		
FD-064(L)-SS	Ø6.4 x 412		
<i>Small femoral implants (short threaded distal fixation)</i>			
FD-032(S)-SS	Ø3.2 x 197		
FD-040(S)-SS	Ø4.0 x 333		
FD-048(S)-SS	Ø4.8 x 404		
FD-056(S)-SS	Ø5.6 x 405		

FEMORAL IMPLANT SPECIFICATIONS

CATALOG #	SIZE Ø x LENGTH [MM]	PROXIMAL FIXATION [MM]	DISTAL FIXATION [MM]
<i>LON femoral implants (non-threaded distal fixation)</i>			
FDLON-F032-SS	Ø3.2 x 194		
FDLON-F040-SS	Ø4.0 x 330		
FDLON-F048-SS	Ø4.8 x 400		
FDLON-F056-SS	Ø5.6 x 411		
FDLON-F064-SS	Ø6.4 x 411	